

### Chamundi Dental Chair Mount Unit

Body Contoured Electrically operated dental chair equipped with state-of-art microprocessor based chair control system.

- Soft key switches for chair operation.
- Two customizable programmes, zeroposition and gargle position.
- Synchronised movement of back rest and seat rest for patient comfort.
- DArticulated Headrest with thickly cushioned rexine.
- · Swivelling arm rest for easy patient access.

## Maia LED (OEM, Faro, Italy) Operating light having multi surface reflector having 30000lux.

- On / Off and intensity control by non-touch sensor with 3 directional movements.
- · Light Fitted on well balanced base mount piller & sturdy arm for operator's convenience.
- Monitor mounting arm

Easily cleanable chair side blue colour oval glass spittoon with auto water connection for spittoon and tumbler, manually moveable towards the patient chest position for convenience

High and Low vacuum motorised non-stop CRS suction.

- Suction power of 350 Ltr / min.
- · Auto drain & Auto flush system.

#### Assistant side Chair controls.

- High & Low vacuum suction cannulas with individual filters collect debris with removable tubing to keep hygiene.
- Beecool LED Light cure unit (Cord less)
- 3-way syringe

## Modular delivery system (Hanging / Over the Patient)

- Airotor points with tubing & nut fitting 2 Nos
- Supreme Micro motor 35,000rpm with digital display of speed
- Supreme optic scaler with five scaling tips
- 3-way syringe
- SS Instrument Tray
- X -Ray viewer.

#### **Multifunction Foot Control**

Joystick function for chair movement, Airotor and Micromotor and scaler Operations

# "COMFORT" Dental operator's stool moving on five castors with chromium plated base.

- Pneumatic piston for up down movement, backrest & seatrest.
- Ergonomically designed hand rest for operator's convenience.



FARO MAIA LED: Even the sunlight affects colour perception: at changing of the seasons or hours of the day we see different colours. Therefore it is

important to have a light source that, at any time, is able to give accurate information to the doctor. All FARO lamps offer intensity adjustment so that the



visual function can capture the tiny details of the operative field with minimal effort, constantly keeping the colour temperature at 5,000 K.

























### CONFIDENT DENTAL EQUIPMENTS PVT. LTD.

Plot No. 17-H, Sector I, Phase II, Bidadi Industrial Area,

Bidadi - 562109, Bangalore (INDIA) GSTIN: 29AAACC9198C1ZF

Ph: +91-80 28023000 | E-mail: info@confidental.org | Web: www.confidental.org